



## Anti-Embolism Stocking System

CAP Stockings are the most advanced anti-embolism stockings available in the world. Combining a unique design, the most precise fitting method, superior manufacturing and quality control, CAP Stockings deliver prescribed graduated compression for the safe and effective prevention of DVT.



Knit-in Femoral Gusset prevents tourniquet effect at upper thigh

Inspection hole on bottom of stocking foot

Non-Constrictive Thigh Band keeps stocking up

Blue Heel™ “floating” design eliminates unwanted pressure on the heel and toe and facilitates easy application

**“Inpatients having surgery should be offered thigh-length graduated compression/anti-embolism stockings from the time of admission to hospital.”**

– Venous Thromboembolism NICE Clinical Guideline 46 2007



# Anti-Embolism Stocking System



## ADDITIONAL FEATURES

- Advanced knitting process insures accurate and controlled compression profile
- Floating heel accommodates varying foot sizes, virtually eliminating unwanted excessive foot and instep compression
- Wicking fibers provide improved patient comfort
- Toe design eliminates toe compression
- Color coded for easy size identification
- Latex free
- Machine washable and dryer safe for convenient patient use after discharge

## Broadest size range available

Carolon stockings are available to accurately fit the broadest patient range in the market with sizes from Small to XXX Large in both knee and thigh lengths. Also available are knee and thigh styles in short, regular and long lengths.

## Easy and accurate fitting

Because compression profiles are based on the ankle compression, fitting procedures are based on an ankle measurement for accurate fit and compression.

Constructed to BS 7672:1993 specifications, CAP Stockings are tested by an independent laboratory for compliance.

## ATS<sup>®</sup>



## Bariatric Adjustable Thigh Stocking

Fits up to a 39" (99cm) thigh circumference

- Adjustable hook and loop design eases application for bariatric patients
- Adjustable thigh allows accurate compression for total hip patients



601 Forum Parkway  
Rural Hall, NC 27045 USA  
www.carolon.com

Toll Free: 800.334.0414  
Fax: 336.969.6999  
Email: carolon@carolon.com

Made in USA  
© 2015 Carolon Company  
00404-CB2