

Health Support[®]

MEDICAL COMPRESSION

- Accurately controlled compression begins at the ankle and reduces in intensity up the leg
- Reduce swelling
- Improve leg health
- Fashionable appearance
- Machine washable – Machine dry
- Satisfaction guaranteed
- Latex free – Made in the U.S.A.

Sheers

- Advanced sheer fabric
- Reinforced toe and floating heel design
- Panty hose with wide comfort band
- Thigh length with silicone non-slip band
- Available in black and beige

Socks

- Advanced soft ribbed fabric
- Reinforced flared toe
- Reciprocated Y-heel for added comfort
- Non-roll top
- Available in black and white

To order, visit
www.CarolonStore.com
or call us toll free at
800.334.0414



Fashionable relief for tired aching legs,
swelling and varicose veins



601 Forum Parkway
Rural Hall, NC 27045 USA
800.334.0414
www.carolon.com





SOCK SHOWN

15-20 mmHg

Light pressure helps to prevent edema and leg discomfort. Also helps to prevent deep vein thrombosis (DVT) in individuals subjected to immobility.

Products Available: Sheer Knee, Sock, Sheer Thigh CT, Sheer Thigh OT, Pantyhose and Pantyhose Maternity



SHEER SHOWN

20-30 mmHg

Moderate pressure helps to prevent more severe swelling and discomfort. Often used after laser varicose vein and venous ablation procedures.

Products Available: Sheer Knee, Sheer Knee with Zipper, Sock, Sock with Zipper, Sheer Thigh OT, Sheer Thigh CT, Pantyhose and Pantyhose Maternity

30-40 mmHg

A high level of pressure used for those with severe vascular insufficiency and for the prevention and treatment of venous ulcers.

Products Available: Sheer Knee, Sheer Thigh OT, 2-Layer Sheer and 2-Layer Sock

40-50 mmHg

Products Available: 2-Layer Sheer and 2-Layer Sock



SHEER SHOWN



S i z i n g C h a r t s



Panty Hose and Thigh Length

1 Ankle	2 Calf	3 Thigh	4 Hip		5 Length	Size	
			Standard	Maternity			
7" - 8" 18 cm - 20 cm	10" - 13" 25 cm - 33 cm	16" - 22" 41 cm - 56 cm	To 40" 102 cm	To 52" 132 cm	< 26" (66cm) > 26" (66cm)	Short Regular	A
8" - 9" 20 cm - 23 cm	12" - 15" 31 cm - 39 cm	19" - 25" 49 cm - 64 cm	To 42" 107 cm	To 56" 142 cm	< 28" (71cm) > 28" (71cm)	Short Regular	B
9" - 10" 23 cm - 25 cm	14" - 17" 36 cm - 44 cm	22" - 29" 56 cm - 74 cm	To 48" 122 cm	To 62" 157 cm	< 30" (76cm) > 30" (76cm)	Short Regular	C
10" - 11" 25 cm - 28 cm	16" - 19" 41 cm - 49 cm	24" - 30" 61 cm - 73 cm	To 51" 130 cm	To 68" 173 cm	< 32" (81cm)	Regular	D
11" - 12" 28 cm - 31 cm	18" - 21" 46 cm - 54 cm	24" - 30" 61 cm - 73 cm	To 55" 140 cm	To 72" 183 cm	< 32" (81cm)	Regular	E



Knee Length Sheer

1 Ankle	2 Calf	3 Length	Size	
7" - 8" 18 cm - 20 cm	10" - 13" 25 cm - 33 cm	< 15" (39 cm) > 15" (39 cm)	Short Regular	A
8" - 9" 20 cm - 23 cm	12" - 15" 31 cm - 39 cm	< 16" (41 cm) > 16" (41 cm)	Short Regular	B
9" - 10" 23 cm - 25 cm	14" - 17" 36 cm - 44 cm	< 17" (44 cm) > 17" (44 cm)	Short Regular	C
10" - 11" 25 cm - 28 cm	16" - 19" 41 cm - 49 cm	< 18" (46 cm) > 18" (46 cm)	Short Regular	D
11" - 12" 28 cm - 31 cm	18" - 21" 46 cm - 54 cm	< 18" (46 cm) > 18" (46 cm)	Short Regular	E
12" - 13" 31 cm - 33 cm	20" - 23" 51 cm - 59 cm	< 18" (46 cm) > 18" (46 cm)	Short Regular	F
13" - 14" 33 cm - 36 cm	22" - 26" 56 cm - 66 cm	< 18" (46 cm) > 18" (46 cm)	Short Regular	G