


CAROLON COUTURE™

Couture compression stockings and socks use advanced high-tech microfiber yarns to provide a soft, silky, luxurious feel on the skin and fashionable for wear with any attire. Couture compression stockings and socks are "shaped to fit", combining optimal stretch fabric with a unique contouring process providing ease of donning and superior wearing comfort.

No squeezing or bunching behind the knees.

Carolon Comfort Knit™ technology utilizes yarns uniquely designed to wick perspiration away from the skin for all-day wearing comfort in any climate.




Floating heel for easy donning.

No fabric ridging at the instep

The CAROLON Couture™ Comfort Advantage

CAROLON Couture™ Comfortable Compression Profile
Smooth Pressure Transition From Ankle To Thigh = Unmatched Wearing Comfort



Stockings fit perfectly and deliver nearly imperceptible yet medically correct compression. In properly graduated stockings, compression is highest at the ankle and gradually reduces up to the top of the stocking.



Antimicrobial - Made with 99.9% pure silver, X-Static® enhances socks with permanent anti-odor and antimicrobial performance. It protects the fabric from the growth of odor-causing bacteria and fungus – keep your feet fresh and odor-free!



Permanent - Will not wash or wear out. Performs for the life of the product.



Regulates Temperature - Natural temperature-regulating properties keep your feet cooler in the summer and warmer in the winter.



Anti-Odor - XT2™ technology uses natural silver to provide the ultimate odor protection for the stocking.



Permanent - Will not wash or wear out. Performs for the life of the product.



See how the magic of silver works:

Watch the X-Static video on your smartphone with your QR Code Reader App, or go to carolon.com and choose the Couture tab to watch the video.

Medical Compression Socks for Men and Women

SIZING CHART

ANKLE	CALF	LENGTH	SIZE
7" - 7 3/4" 18 cm - 20 cm	10" - 13" 25 cm - 33 cm	< 15" (39 cm)	A Short
		> 15" (39 cm)	A Regular
8" - 8 3/4" 20 cm - 22 cm	12" - 15" 31 cm - 39 cm	< 16" (41 cm)	B Short
		> 16" (41 cm)	B Regular
9" - 9 3/4" 23 cm - 25 cm	14" - 17" 36 cm - 43 cm	< 17" (43 cm)	C Short
		> 17" (43 cm)	C Regular
10" - 10 3/4" 25 cm - 27 cm	16" - 19" 41 cm - 49 cm	< 18" (46 cm)	D Short
		> 18" (46 cm)	D Regular
11" - 11 3/4" 28 cm - 30 cm	18" - 21" 46 cm - 54 cm	< 18" (46 cm)	E Short
		> 18" (46 cm)	E Regular

CUSHION FOOT SOCKS

Available in Black, Khaki, White

DRESS SOCKS

Available in Black, Brown, Khaki, Navy

COMPRESSIONS

15-20 mmHg, 20-30 mmHg, 30-40 mmHg

Medical Compression Stockings for Men and Women

SIZING CHART

ANKLE	CALF	THIGH	HIP	WAIST	KNEE LENGTH	THIGH LENGTH	SIZE
7 1/2" - 8 1/4" 19 cm - 21 cm	11" - 15" 28 cm - 38 cm	18" - 23" 45.5 cm - 58.5 cm	35" - 42" 89 cm - 107 cm	25" - 32" 63.5 cm - 81 cm	< 16" (41 cm)	< 28" (71 cm)	A Short
					> 16" (41 cm)	> 28" (71 cm)	A Regular
8 1/2" - 9 3/4" 21.5 cm - 25 cm	13" - 17" 33 cm - 43 cm	21" - 26" 53.5 cm - 66 cm	41" - 48" 104 cm - 122 cm	31" - 38" 78.5 cm - 96.5 cm	< 17" (43 cm)	< 29" (74 cm)	B Short
					> 17" (43 cm)	> 29" (74 cm)	B Regular
10" - 11 1/4" 25.5 cm - 28.5 cm	15" - 19" 38 cm - 48.5 cm	24" - 29" 61 cm - 73.5 cm	47" - 54" 119.5 cm - 137 cm	37" - 44" 94 cm - 112 cm	< 17" (43 cm)	< 29" (74 cm)	C Short
					> 17" (43 cm)	> 29" (74 cm)	C Regular
11 1/2" - 13 1/2" 29 cm - 34 cm	17" - 21" 43 cm - 53 cm	27" - 32" 68.5 cm - 81 cm	53" - 60" 134.5 cm - 152 cm	43" - 50" 109 cm - 127 cm	< 17" (43 cm)	< 29" (74 cm)	D Short
					> 17" (43 cm)	> 29" (74 cm)	D Regular

STOCKINGS

Available in Black, Beige

COMPRESSIONS

20-30 mmHg, 30-40 mmHg

STYLES

Below Knee, Full Length Thigh with Silicone Dot Stay Up Band, Full Length Thigh with Waist Attachment, Tights, Maternity Tights

CARE INSTRUCTIONS

Couture stockings and socks should be laundered daily for maximum effectiveness and to prolong the life of the garments.

MACHINE WASHING (Recommended)

Machine wash in hot water. Dryer safe on medium heat. Do not use softeners.

Periodically wiping the silicone dots inside the top band with an alcohol swab will remove the accumulation of natural body oils and lotions ensuring the stocking will remain securely in place throughout the day.



Do not use chlorine based bleach.
Do not iron.
Do not dry clean.

THE COUTURE TOTAL SATISFACTION GUARANTEE

If, for any reason, you are not pleased with Couture, simply wash the product and return it for a replacement or refund.



To place an order or for more information, please call 877.232.8825 or fax Customer Service 877.232.8855.



601 Forum Parkway
Rural Hall, NC 27045 USA
www.carolon.com

European Authorized Representative:
Emergo Europe
Molenstraat 15 • 2513 BH
The Hague, Netherlands