



Multi-Layer Stocking System

Advanced Compression System for the Treatment of Venous Leg Ulcers



Precise Compression

Overstocking and Understocking are designed to work together to provide clinically appropriate compression. Available in 30-40mmHg (A6531) or 40-50mmHg (A6532). Understocking can be worn alone during sleep to provide 15-20mmHg while patient is recumbent.

Ease of Application

Application and removal are more easily facilitated since compression is divided among two layers. The fabric on the Overstocking has been specifically designed to slide comfortably over the Understocking.

A Transition to Healing

After the patient's edema has been stabilized with the Carolon Wrap System, the Multi-Layer Stocking system should be used.

The advanced Multi-Layer Stocking provides accurate compression, proven healing performance, superior patient compliance and allows the patient to take part in the application of their own compression therapy.

3-Step Carolon Compression System

1



For the initial treatment, apply the Carolon Wrap System until the patient is stabilized.

2



For the final phase of compression treatment, use the Carolon Multi-Layer Hosiery System.

3



To help prevent recurrence, Carolon has a full line of medical compression hosiery.



601 Forum Parkway
Rural Hall, NC 27045 USA
www.carolon.com

European Distributor: H&R Healthcare Ltd.
3 Redcliff Rd., Melton Park
Melton, Hull UK HU14 3RS

European Authorized Representative:
Emergo Europe
Molenstraat 15 • 2513 BH
The Hague, Netherlands


800.334.0414
Made in USA

© 2014 Carolon Company



Multi-Layer Stocking System

Smallest Ankle Circumference	Largest Calf Circumference	Beige (Sheer Fabric)	Black (Sock Fabric)	Length	Size	Beige	Black
		Bottom of heel to bend in knee					
18cm - 20cm	25.5cm - 33cm	To 38cm	To 40.5cm	Short	A	81001122	82601042
		Over 38cm	Over 40.5cm	Regular		81011122	82611042
20cm - 22.5cm	30.5cm - 38cm	To 40.5cm	To 43cm	Short	B	81002122	82602042
		Over 40.5cm	Over 43cm	Regular		81012122	82612042
22.5cm - 25.5cm	35.5cm - 43cm	To 43cm	To 45.5cm	Short	C	81003122	82603042
		Over 43cm	Over 45.5cm	Regular		81013122	82613042
25.5cm - 28cm	40.5cm - 48cm	To 45.5cm		Short	D	81004122	82604042
		Over 45.5cm		Regular		81014122	82614042
28.5cm - 30.5cm	45.5cm - 53cm	To 45.5cm		Short	E	81005122	82605042
		Over 45.5cm		Regular		81015122	82615042
30.5cm - 33cm	50cm - 58.5cm	To 45.5cm		Short	F	81006122	N/A
		Over 45.5cm		Regular		81016122	N/A
33cm - 35.5cm	55.5cm - 66cm	To 45.5cm		Short	G	81007122	N/A
		Over 45.5cm		Regular		81017122	N/A



Measurements:

1. Smallest ankle circumference
2. Largest calf circumference
3. Length from base of heel to bend in knee
4. Use chart above to determine size and order number

To order:

Please contact your DME provider.



601 Forum Parkway
Rural Hall, NC 27045 USA
www.carolon.com

European Distributor: H&R Healthcare Ltd.
3 Redcliff Rd., Melton Park
Melton, Hull UK HU14 3RS

European Authorized Representative:
Emergo Europe
Molenstraat 15 • 2513 BH
The Hague, Netherlands



800.334.0414

Made in USA

Rev: 00404-ML1